

Directionality Characteristics of the Tympan Open-Source Hearing Aid and Earpieces

Benjamin P. Masters

Joshua M. Alexander


PURDUE
UNIVERSITY

Methods

- When possible, followed procedures outlined in **ANSI/ASA S3.35-2021**
 - *“Method of measurement of performance characteristics of hearing aids under simulated real-ear working conditions”*
- **Pink Noise Stimulus**
 - ≈ 85 dB SPL at hearing aid mic location
 - 5-second duration
 - Analyzed in 1/3-octave bands


Set up

- KEMAR centered in a double-walled sound booth (2.13 m x 2.44 m)
 - Rotated in 10° increments
 - **Left earpiece was tested**
- Hafler M5 Reference Monitor sound field loudspeaker
 - Flat frequency response




Earpieces

- 9.5-mm mic spacing (d)
 - $c \approx 343$ m/s
 - Travel time (d/c): $27.7 \mu\text{s}$
- Sampling rate set to 36.1 kHz \rightarrow 1-sample delay = $d/c = 27.7 \mu\text{s}$
 - Cardioid directivity



Microphone Matching

- Critical for optimal cancellation
 - **Rear mic output** was ≈ 5 dB lower than the **front mic**
 - 5 dB of gain was added to **match the two mics**



Omni-Directional Responses

- It is well-known that putting mics on the head will **shift the response** to the side due to head shadow


Front Mic (re 0°)

Tympan Directionality - Left Side

Front Mic

0 dB Gain

AIDI = -1.12




Rear Mic (re 180°)

Tympan Directionality - Left Side

Rear Mic

0 dB Gain

AIDI = -1.00




Rear Mic Delay

→
0° reference

-3 dB Rear Gain

Tympan Directionality - Left Side
1 sample Rear Delay
-3dB Rear Gain


AIDI = +2.31
(+3.43)



0 dB Rear Gain

Tympan Directionality - Left Side
1 sample Rear Delay
0dB Rear Gain


AIDI = +0.50
(+1.62)



+3 dB Rear Gain

Tympan Directionality - Left Side
1 sample Rear Delay
+3dB Rear Gain

AIDI = -1.53
(-0.41)



*Gain values after mic matching

AIDI was lowest with -6 dB and +6 dB rear mic gain

AIDI = Articulation Index weighted Directivity Index
values in () = Improvement re: front mic only


Front Mic Delay

←
180° reference

-3 dB Front Gain

Tympan Directivity - Left Side
1 sample Front Delay
-3dB Front Gain


AIDI = +0.95
(+1.95)



0 dB Front Gain

Tympan Directivity - Left Side
1 sample Front Delay
0dB Front Gain


AIDI = +2.98
(+3.98)



+3 dB Front Gain

Tympan Directivity - Left Side
1 sample Front Delay
+3dB Front Gain

AIDI = -0.11
(+1.11)




*Gain values after mic matching

AIDI was lowest with -6 dB and +6 dB front mic gain

AIDI = Articulation Index weighted Directivity Index
values in () = Improvement re: rear mic only


AIDI Improvement

- Front mic with rear mic delay and gain = -3 compared to front mic only (0° reference)
- Rear mic with front mic delay and gain = 0 compared to rear mic only (180° reference)




Future Directions: Variable Directionality

- **Step 1:** create opposing cardioid patterns
 - Front facing
 - Rear facing: inverted phase and variable gain, b





(Kates 2008)

Future Directions: Variable Directionality




The **front** and **rear** responses **add with opposite sign**, thereby creating nulls where they intersect

- **Step 2:** vary relative gain b on rear cardioid to create different polar patterns with a continuous range of nulls



$b = 0.5$



$b = 0.25$

Conclusions

1. Optimized directionality requires **precise timing** (sampling rate) and **microphone matching**
2. Directional patterns and relative AIDI improvements ($\approx 3\text{-}4$ dB) are **consistent with published research**
3. **Optimal directivity**
 - Front cardioid obtained with -3 dB gain on the rear mic
 - Rear cardioid obtained with matched mic gain
4. **Future directions**: use optimized front and rear cardioids with opposite phase to **vary null with fixed and adaptive directionality**