



# Immersive Multitalker Remote Microphone System

Ryan M. Corey and Andrew C. Singer

181<sup>st</sup> Acoustical Society of America Meeting  
Seattle, Washington      December 2, 2021



UNIVERSITY OF  
**ILLINOIS**  
URBANA-CHAMPAIGN

# Remote Microphones



**Problem:** Hearing aids work poorly in noise

Remote microphone accessories:

- Transmit sound directly from talker to listener
- Best way to improve intelligibility in noise!

Limitations of current systems:

- Limited to one talker at a time
- Mono signal – no spatial cues or room acoustics




## Goal: Immersive listening system for multiple talkers

- Multiple simultaneous talkers
- Preserve spatial cues
- Low background noise
- Track motion

### Key idea:

Combine low noise of remote mics with spatial cues of earpiece mics



### Remote microphones

- ✗ No spatial cues
- ✓ Low noise

### Earpiece microphones

- ✓ Spatial cues
- ✗ Noisy


### Proposed system

- ✓ Spatial cues
- ✓ Low noise

# Adaptive Binaural Filter



**Solution:** Filter remote mic signals to match magnitude & phase at earpiece mics



Filter adaptively minimizes error signal  $|y - x_{\text{ear}}|^2$

R. Corey and A. Singer, "Adaptive binaural filtering for a multiple-talker listening system using remote and on-ear microphones," *IEEE Workshop on Applications of Signal Processing to Audio and Acoustics (WASPAA)*, October 2021.



# Tympan Implementation



## 2x Wireless mic transmitters

Saramonic UwMic9

Near-zero latency with analog link!

## Tympan with AIC CODEC Shield

4 in / 4 out

## Earbuds with binaural mics

Roland CS-10EM


Stereo in / stereo out

## Stereo wireless mic receiver


Mic A → “left” channel

Mic B → “right” channel

# Tympan Implementation



- Rechargeable battery
  - First fully portable prototype I've built!
- Android app
  - User control of mixer levels
  - Performance monitoring
- Multichannel recording to SD card
  - Compare processed & unprocessed binaural signals
  - Quantify filter performance using known signals





Unprocessed  
Processed




## Binaural demo

Best experienced  
with headphones

[go.illinois.edu/  
AugmentedListening](https://go.illinois.edu/AugmentedListening)



- Adaptive filtering improvements
  - Parameterized model using fewer coefficients
  - Refine adaptive filter for faster motion tracking
  - Multiple-input filter for arrays (see WASPAA paper)
- Hearing aid features
  - Frequency-selective gain, feedback control
  - Per-talker dynamic range compression
- Tympan BTE earpieces
  - Modify to work with external mics



[github.com/Tympan](https://github.com/Tympan)



[go.illinois.edu/AugmentedListening](https://go.illinois.edu/AugmentedListening)



@ryanmcorey

Supported by:

- Intelligence Community Postdoctoral Research Fellowship Program
- National Science Foundation Partnership for Innovation Program
- Discovery Partners Institute Research and Development Program